



# HYSYNAL™ SYNTHETIC RUBBER TOURNIQUETS

Hygenic's knowledge and expertise in elastomers played an integral part in introducing one of the first widely available, non-latex disposable tourniquets to the healthcare industry. Disposable tourniquets have become a staple in the healthcare industry by helping to prevent the spread of infection and contamination. For the last 35 years, clinics, hospitals and healthcare professionals have chosen HYSYNAL Synthetic Rubber Tourniquets as their preferred tourniquet strap.

## PRODUCT OVERVIEW

- Proprietary technology creates non-slip texture that ensures the tourniquet stays tight and tied, resulting in more efficient needle sticks
- Exclusive polymer formulation eliminates concerns over latex proteins
- Our synthetic polyisoprene compound provides better tear resistance and faster recovery, as compared to TPE
- Powder free
- Made in the USA
- Standard color is blue
- Standard sizes are ¾" or 1" wide by 18" long

## LOOKING FOR A UNIQUE, CUSTOM TOURNIQUET?

Hygenic can also create custom tourniquet straps to suit your specific needs. We have the ability to produce custom colors and sizes. We can also print logos, symbols or words on the tourniquet straps.

## AVAILABLE PACKAGING OPTIONS

Hygenic has a variety of standard packaging options that suit most customers' needs. Standard tourniquets are sold in quantity batches of 1,000 straps.

### 1" TOURNIQUETS

#### LAY FLAT

- 10 boxes of 100 straps
- 4 boxes of 250 straps

#### ON A REEL

- 25 perforated straps/reel with 40 reels/case


### ¾" TOURNIQUETS

#### BULK

- 1,000 straps dropped into a box

## NEED A DIFFERENT PACKAGING SOLUTION?

We can also create custom packaging solutions for special order.



For more details or to place an order, contact us online or by phone.

[www.hygenic.com/contact-us](http://www.hygenic.com/contact-us)  
1-800-321-2135